



高新技术企业

Hi-Tech Enterprise




美國榮佳醫療有限公司

WickiMed Supply Company Ltd.


Welcome 欢迎

ようこそいらっしゃいませ


Group Structure



Milestone



Milestone



Equipment

Injection Molding Workshop – 300K cleanroom


6 sets
Vertical
Injection
Machine

6 sets
Horizontal
Injection
Machine



Injection Molding Components

Double Injection



PC / ABS / PE / PP



Equipment

Turning / Stamping / Machining

20 CNC Turning /
Machining

55 Stamping
Machine

40 Turning
Machine



Equipment Tooling Marking

8 Grinding Machine



2 EDM Machine



3 Wire Cut Machine

Countries & Markets Served

Domestic

93 cities in Mainland China

Overseas

USA & Europe

- * Italy, Austria, Spain, France
- * USA, Germany, France
- * UK, Ireland

Asia

- * Israel, Iran, India
- * Korea, Hong Kong
- * Taiwan, Malaysia
- * Philippines, Japan
- * Australia


Products - Electrosurgical Pencil / Suction Coagulator

Electrosurgical Pencils


Hand-switch / foot control electrosurgical pencil

Main feature

- Suitable for different kinds of surgical
- Closed pencil body
- Hexagon locking system
- Changeable electrode tips
- Non-slip pencil body




P10S
Hand-switch electrosurgical pencil
Fitting with standard blade tips, sterilization packaging



P10S
Hand-switch electrosurgical pencil
Fitting with standard blade tips, sterilization packaging



F18S
Foot control electrosurgical pencil
Fitting with standard blade tips, sterilization packaging



P10S
Hand-switch electrosurgical pencil
Fitting with standard blade tips, sterilization packaging, 1pcs/bag, 10pcs/case, 100pcs/carton



P10SC
Hand-switch electrosurgical pencil
Fitting with standard blade tips, sterilization packaging, 1pcs/bag, 10pcs/case, 100pcs/carton



P11SC
Hand-switch electrosurgical pencil
Fitting with needle tips, sterilization packaging, 1pcs/bag, 10pcs/case, 100pcs/carton



P30S
Hand-switch electrosurgical pencil
Fitting with standard blade tips, sterilization packaging, 1pcs/bag, 5pcs/case, 50pcs/carton

P10R
Hand-switch electrosurgical pencil
Fitting with standard blade tips, sterilization packaging, 1pcs/bag, 6pcs/case, 30pcs/carton

F18S
Foot control electrosurgical pencil
Fitting with standard blade tips, sterilization packaging, 1pcs/bag, 10pcs/case, 100pcs/carton

P10SX
Hand-switch / foot control suction coagulator electrosurgical pencil
Fitting with extended blade tip, sterilization individual packaging

P20TS
Hand-switch electrosurgical pencil with holster
Sterilization individual packaging, 60set/carton

Provide three-in-one combination packaging, to meet your need
Electrosurgical pencil + grounding pad + tip cleaner



Products - Electrosurgical Pencil

Electrosurgical Pencils |

Telescopic / rotational telescopic electrosurgical pencil

Main features

- Telescopic electrode tip deeper into the tissue without the need to change the head which is more flexible and efficient
- Flexible distance, distance between electric switch and electrode head the longest extended distance is 130mm
- The humanity design, perfect and comfortable line non-slip body is more comfortable and convenience
- Suitable for the department of cardiac surgery, gynecology, gastrointestinal surgery, hepatopancreatobiliary surgery, laparoscope and orthopedics



Y90S

Telescope electrosurgical pencil
with blade electrode tip,
sterilization packaging

Y91S

Telescope electrosurgical pencil
with needle electrode tip,
sterilization packaging

Y92S

Telescope electrosurgical pencil
with L hook electrode tip,
sterilization packaging

Y93S

Telescope electrosurgical pencil
with tungsten needle electrode tip,
sterilization packaging

WickiMed®

Y090K

Rotational telescopic electrosurgical pencil
with blade electrode tip,
sterilization packaging

Y091K

Rotational telescopic electrosurgical pencil
with needle electrode tip,
sterilization packaging

Y092K

Rotational telescopic electrosurgical pencil
with L hook electrode tip,
sterilization packaging

Y093K

Rotational telescopic electrosurgical pencil
with tungsten needle electrode tip,
sterilization packaging


Products - Electrosurgical Pencil

Electrosurgical Pencils

- Electrosurgical pencil with L-hook / Suction electrosurgical pencil with L-hook

Main features

- Suitable for thoracic surgery, general surgery, gynecology, minimally invasive surgery and other endoscopic surgery
- Effectively to separate tissue, coagulation and cutting, is stable and safety
- Non-stick coated electrode tip, can reduce tissue adhesion during cutting surgery
- Smoke valve, can evacuate smoke, to protect the medical workers' health



T330K
Electrosurgical pencil
with L-hook,
sterilization packaging

T350K
Electrosurgical pencil
with L-hook,
sterilization packaging

Y330K
Smoke evacuation
electrosurgical pencil
with L-hook,
sterilization packaging


Y350K
Smoke evacuation
electrosurgical pencil
with L-hook,
sterilization packaging

Electrosurgical Pencils

- LED electrosurgical pencil

Main features

- Suitable for different kinds of surgery
- Adjust the length as required
- Closed body design
- Non-slip body design
- With LED light



D010K
LED electrosurgical pencil
with blade electrode tip,
sterilization packaging

D011K
LED electrosurgical pencil
with needle electrode tip,
sterilization packaging

D100S
LED electrosurgical pencil
with blade electrode tip,
sterilization packaging


Products - Smoke Pencil

Electrosurgical Pencils |

● Multifunctional telescopic smoke evacuation electrosurgical pencil

Main features

- Hand switch smoke evacuation function
 - The smoke evacuation tube and telescopic electrode tip with locked device
 - The electrode tip can be changed to any types and adjust length as required by the operation, the longest expended distance is 35cm
 - Suitable for the department of cardiac surgery, gynecology, gastrointestinal Surgery, hepatopancreatobiliary surgery clinic, laparoscope, orthopedics and etc.
- A** Telescopic property is suitable for surgery required deep operation
Adjust the length at any time according to the depth of the surgical field of vision, greatly enhancing the flexibility of clinical use
- B** The smoke evacuation tube can directly and quickly suck away the smoke and exhaust gas during the operation
- C** Can directly connect with smoke ejector for smoke evacuation
The transparent medical PC air channel is thread design, it is stably and conveniently to change the electrode tip



- | | | | | |
|---|--|--|---|--|
| <p>Y70S
Telescopic smoke evacuation electrosurgical pencil with blade electrode tip, sterilization packaging</p> | <p>Y71S
Telescopic smoke evacuation electrosurgical pencil with needle electrode tip, sterilization packaging</p> | <p>Y72S
Telescopic smoke evacuation electrosurgical pencil with L hook electrode tip, sterilization packaging</p> | <p>Y73S
Telescopic smoke evacuation electrosurgical pencil with tungsten needle electrode tip, sterilization packaging</p> | <p>AF4797SE
Smoke evacuation electrosurgical pencil with blade electrode tip, sterilization packaging</p> |
|---|--|--|---|--|

Products - Grounding Pad



WG300-R20
With REM plug

WG300-R30
without REM plug

WG300-R40
with HFI 6.3 plug

Adult/Pediatric grounding pad

WG300-C20

Disposable bipolar adult/pediatric GP without pre-wire

WG300-C25

Disposable monopole adult/pediatric GP without pre-wire

Pediatric grounding pad

WG300-P10

Disposable bipolar pediatric GP without pre-wire

WG300-P15

Disposable monopole pediatric GP without pre-wire

WG300-P20

Disposable bipolar pediatric GP with REM plug

WG300-P30

Disposable monopole pediatric GP without REM plug

Adult grounding pad

WG300-A10

Disposable bipolar adult GP without pre-wire

WG300-A15

Disposable monopole adult GP without pre-wire

WG300-A20

Disposable bipolar adult GP with REM plug

WG300-A30

Disposable monopole adult GP with HFI 6.3 plug



Grounding Pad

Main features

- Suitable for different kinds of surgery
- Special anti - allergy adhesive layer proof firm attachment, safe and comfortable
- Provide adult, pediatric and general series



Products - Series of Electrode



Electrode
Blade / Ball / Needle

Blade Electrode |



Ball Electrode |



Needle Electrode |



Blade electrode	Ball electrode	Needle electrode
WP200-B100 Blade	WP200-C100 Ball 3mm	WP200-N100 Needle 0.65mm
WP200-B101T Blade (with heat shrinkable tube)	WP200-C102 Ball 4mm	WP200-N101 Needle 0.80mm
WP200-B203T Extended blade(10cm /Teflon)	WP200-C101 Ball 5mm	WP200-N300 Angle needle 0.80mm
WP200-B300 Angle blade	WP200-C300 Ball 3mm/45°Angle	WP200-N200 Extended needle 0.65mm
WP200-B200 Extended blade	WP200-C301 Ball 5mm/45°Angle	WP200-N201 Extended needle 0.80mm
WP200-B200T Extended blade(Teflon)	WP200-C200 Extended ball 3mm	WP200-N301 Extended angle needle 0.80mm
WP200-B202 Extended blade (with heat shrinkable tube)	WP200-C202 Extended ball 4mm	WP200-N100RN Reusable needle
WP200-B102RN Reusable blade	WP200-C201 Extended ball 5mm	WP200-N102 Straight tungsten needle 0.6mm
WP200-B200RN Reusable extended blade	WP200-C203RN Reusable extended ball 4mm	WP200-N302 Angle tungsten needle 0.8mm



WICKIMED SUPPLY COMPANY LTD

Products - Series of Electrode

LAPAROSCOPIC ELECTRODE

● NON-STICK COATING TECHNOLOGY

REDUCE THE SCAB | LESS CLEANING

DESIGNED FOR PROFESSIONAL SURGICAL NEEDS

OPTIMAL PERFORMANCE OF LAPAROSCOPIC SURGERY

SUMMARIZE |

The laparoscopic electrode introduces non-stick coating technique, with the advantages of a on-stick coating compared to the non-coating:

- Reduce Sticky
- Reduce Eschar

ADVANTAGE |

The electrode polytetrafluoroethylene (PTFE) provides the effect of a non-stick coating without affecting operational performance.

The reduction of eschar accumulation can greatly reduce the frequency of electrode entering and leaving the casing, and improve the surgical efficiency.

The laparoscopic electrode makes the fine operation in the high surgical process, which can effectively reduce the accumulation of eschar and cause little damage to human tissue.

EFFECT |

Apply with non-adhesive material (PTFE) for cutting and fixing soft tissue. The PTFE coating reduces the accumulation of eschar on the electrode during use; therefore, eliminates the need to remove eschar accumulation on the scraping abrasive surface during cleaning.

SPECIFICATIONS |

Working length up to 45CM, including various head shapes (Standard Blade, Round Spatula, Curved Blade, L-hook, J-hook, L-wire, etc.)

WK WickiMed®



LAPAROSCOPIC ELECTRODE

● NON-STICK COATING TECHNOLOGY

REDUCE THE SCAB | LESS CLEANING

LAPAROSCOPIC CURVED BLADE ELECTRODE




WK-U0019
WK-U0019L
WK-U0019S

LAPAROSCOPIC L-HOOK ELECTRODE



WK-U0020
WK-U0020L
WK-U0020S



355mm

475mm

WK-U0019

WK-U0019S

WK-U0020

WK-U0020S

WK-U0019L

WK-U0020L

Products - Series of Electrode

LAPAROSCOPIC ELECTRODE

● NON-STICK COATING TECHNOLOGY

REDUCE THE SCAB | LESS CLEANING

LAPAROSCOPIC J-HOOK ELECTRODE




WK-U0021
WK-U0021L
WK-U0021S

LAPAROSCOPIC L-WIRE ELECTRODE



WK-U0100
WK-U0100L
WK-U0100S



WK-U0021
WK-U0021S
WK-U0100
WK-U0100S
WK-U0021L
WK-U0100L

LAPAROSCOPIC ELECTRODE

● NON-STICK COATING TECHNOLOGY

REDUCE THE SCAB | LESS CLEANING

LAPAROSCOPIC STANDARD BLADE ELECTRODE




WK-U0017
WK-U0017L
WK-U0017S

LAPAROSCOPIC ROUND SPATULA ELECTRODE



WK-U0018
WK-U0018L
WK-U0018S




WK-U0017
WK-U0017S
WK-U0018
WK-U0018S
WK-U0017L
WK-U0018L

Products - LEEP Tungsten Electrode

LEEP Electrode

Main feature

- It is wide application range of surgery operation, especially suitable for the minimal surgery with high precision requirement.
- Various shapes of electrodes provide the best effect in cutting and homeostasis.
- Minimize cell destruction; it gives good cosmetic result and shortens the operation time and healing period.
- It provides wide range of functions including precutting, coagulation, hemostatic, skin incision and biopsy.



WickiMed®



WICKIMED SUPPLY COMPANY LTD

LEEP electrode

€ 2.36



E1015
LEEP needle



E3015
LEEP angle needle



E2015
LEEP ball electrode € 2.3mm



E955
LEEP rhombus electrode



E965
LEEP triangle tungsten electrode



E1025
LEEP round electrode € 4.00



E945
LEEP round electrode € 6.65



E1035
LEEP round electrode € 8.00



E1045
LEEP round electrode € 10.00

LEEP electrode

€ 1.58



E905
LEEP ball needle



E915
LEEP rhombus electrode



E925
LEEP triangle electrode



E935
LEEP angle needle electrode



E1005
LEEP needle electrode



E2005
LEEP ball electrode



E3005
LEEP bending needle electrode

WickiMed®


Products - Loop Tungsten Electrode

Tungsten Electrode

Loop

Main features

- Fitted with all kinds of electrosurgical pencil
 - Precision 'CUT' effect and excellent 'COAG' function
 - Reduce the recurrence rate of pathological tissue
 - Stable 'Y' construction design secured the fine wire Loop
 - Tungsten wire suits continue cutting operation
- More colors and styles can be customized ,OEM



Products - Loop Tungsten Electrode

Tungsten Electrode

Square / Triangle / Halberd / B-Shape

Square tungsten electrode



Main Features:

Fitted with all kinds of electrosurgical pencil
Precision 'CUT' effect and excellent 'COAG' function
Reduce the recurrence rate of pathological tissue
Stable 'Y' construction design secured the fine wire Loop
Tungsten wire suits continue cutting operation

Triangle/ halberd/ B-shape tungsten electrode

Main Features:

Fitted with all kinds of electrosurgical pencil
Precision 'CUT' effect and excellent 'COAG' function
Reduce the recurrence rate of pathological tissue
Stable 'Y' construction design secured the fine wire Loop
Tungsten wire suits continue cutting operation



Square Tungsten Electrode

Triangle / Halberd / B-Shape Tungsten Electrode



Products - Bio-Cone (Patented)

BIO-CONE

Rotational Cone Biopsy Electrode





Easy To Use And Safety For Patients

- A very easy method available to perform cone biopsy in an office or clinic under local or general anesthesia with less or no bleeding
- Controlled movement of wire produces a clean and smooth cut with a minimum of thermal influence. The tissue is conveniently secured for pathological evaluation, up to the cutting edges
- Insulated components ensure safety to the patient, physician and staff
- Stabilizing prongs in combination with rotational hub make sure satisfaction in comfort and safety
- Precision 'CUT' effect to avoid the risk of removing the normal tissue
- Reduce the recurrence rate of pathological tissue
- Fitted with all kinds of electrosurgical pencil



Conization At The Portion Uteri With The BIO-CONE Electrode



After cone biopsy with the BIO-CONE electrode the patients feel better with less pain, less bleeding!

1. Activate the electrosurgical unit; insert the appropriately sized BIO-CONE electrode into the endocervical canal
2. After the stabilizing prongs have secured the cervix, gently rotate the shaft spindle for one complete rotation. The hub serves as rotation guide.
3. Deactivate the electrosurgical unit now and remove the BIO-CONE electrode with the secured sample.


Products - Tungsten Needle

Micro Tungsten Needle

Main features

- Precise tissue dissection and control bleeding, speedy healing through bloodless operation
- Output current kept at minimal as the micro contacting point is at the needle tip.
- Precise cutting and micro-carving to minimize tissue adhesion and prevent bleeding
- Neat and tidiness operation thus shorten the surgical duration

Straight needle




long ace-shaped needle



WickiMed®



WICKIMED SUPPLY COMPANY LTD



WickiMed®


Products - Forceps

Reusable Bipolar Forceps

Main feature

- The tip of forceps processing with special material, with excellent anti- viscosity
- Provide different width of forceps tip for different surgery
- Bayonet and straight design is suitable for different specialized surgery
- Provide different forceps length as surgery required.
- Handle elastic moderate, feel good, long time operation is not easy to fatigue

Wk WickiMed®



Handle elastic moderate


Long time operation is not easy to fatigue

The tip of forceps processing with special material

With excellent anti- viscosity

FDA CFDA CE 1639 PMDA  

WICKIMED SUPPLY COMPANY LTD



BF-01S05110

Straight

Straight copper alloys tip, length 110mm

BF-01S05130

Straight

Straight copper alloys tip, length 130mm

BF-01S05160

Bayonet

Straight copper alloys tip, length 160mm

BF-01B05185

Bayonet

Straight copper alloys tip, length 185mm

BF-01B05200

Bayonet

Straight copper alloys tip, length 200mm

BF-01B10220

Bayonet

Straight copper alloys tip, length 220mm

BF-07B05200

Bayonet

Straight stainless steel tip, length 200mm

BF-07B05210

Bayonet

Straight stainless steel tip, length 210mm

BF-02S05130

Straight

Straight stainless steel tip, length 130mm

BF-05B10210

Bayonet

Straight stainless steel tip, length 210mm

BF-06B05210

Bayonet

Angle stainless steel tip, length 210mm

Wk

35/36
WICKIMED SUPPLY

Products - Veress Needle

7 WickiMed Supply Company LTD

Trocar

Comparison of specification and model



► Specifications	Cannula DIA (mm)	Cannula length(mm)
Bladeless Trocar	5/10/12/15	75/100/150
Bladed Trocar	5/10/12/15	75/100/150
Optical Trocar	5/10/12/15	75/100/150

8

Veress Needle

- This needle is composed of two stainless steel tubes. The outer tube is sharp for penetrate and the inside tube is for air inflation. Also, there is an indicator in the inside tube. When the outer tube penetrate the cortex, the red indicator will move down, which means the penetrating has been completed. Simple and easy to use!
- This needle is packed with non-toxic, aseptic and disposable package, which reduce the risk of infection by crossover operation.



Self-lock structure



Indicator Device



Products - Bladeless Trocar

1 WickiMed Supply Company LTD

2

Bladeless Trocar



Products features:

Before and after the laparoscopic minimally invasive surgery, the puncture guide tube will create a barrier-free access in abdomen so as to leave enough space for other equipments. The abdomen needs keeping inflated and maintaining pressure, in order that there is enough operation space between viscera and skin and pretend incidental injury to viscera.

■ Excellent dynamic and static tightness



A To avoid pressure leakage, the round-hole silicon sheet can attach the equipments during the processing of abdominal puncture or other

B To avoid pressure leakage, the funnel-shaped slotted-cut silicon sheet can be combined automatically after the trocar and other equipments pull out from the guide tube.

■ Strong fixity of abdomen



The cannula is thread design and inverted triangle shape, in order to increase the pull-out resistance and insure that the external wall of guide tube fit tightly with the abdominal wall and not easy off.

■ Combination of separate parts



To solve the problem of keeping the abdomen inflated and pressure maintenance during the operation, the puncture guide tube will adopt the way of combining the round-hole silicon sheet and the funnel-shaped slotted-cut silicon sheet, which can respectively bear 4Kpa pressure.

Products - Bladed Trocar

3 WickiMed Supply Company LTD

Bladed Trocar



Concealed knife metal head



Under the operation of key press and hand press, the concealed knife metal head will act on the surgical site and help smooth the puncture. The keylock mode is able to prevent accidental operation, which efficiently ensures the safety of doctors and patients.

4

Excellent dynamic and static tightness



- A To avoid pressure leakage, the round-hole silicon sheet can attach the equipments during the processing of abdominal puncture or other
- B To avoid pressure leakage, the funnel-shaped slotted-cut silicon sheet can be combined automatically after the trocar and other equipments pull out from the guide tube.

Strong fixity of abdomen



The cannula is thread design and inverted triangle shape, in order to increase the pull-out resistance and insure that the external wall of guide tube fit tightly with the abdominal wall and not easy off.

Combination of separate parts



To solve the problem of keeping the abdomen inflated and pressure maintenance during the operation, the puncture guide tube will adopt the way of combining the round-hole silicon sheet and the funnel-shaped slotted-cut silicon sheet, which can respectively bear 4Kpa pressure.

Products - Optical Trocar

5 WickiMed Supply Company LTD

6

Optical Trocar



Optical Trocar



In the complex abdominal wall and abdominal environment, the visual trocar can be visually inserted into the abdominal cavity through abdominal wall with the help of camera. In this case, doctors can detect puncture accidents in time and reduce the risk of wounding abdominal wall tissue and abdomen organs by accidental operation.

Excellent dynamic and static tightness



- A To avoid pressure leakage, the round-hole silicon sheet can attach the equipments during the processing of abdominal puncture or other
- B To avoid pressure leakage, the funnel-shaped slotted-cut silicon sheet can be combined automatically after the trocar and other equipments pull out from the guide tube.

Strong fixity of abdomen



The cannula is thread design and inverted triangle shape, in order to increase the pull-out resistance and insure that the external wall of guide tube fit tightly with the abdominal wall and not easy off.

Combination of separate parts



To solve the problem of keeping the abdomen inflated and pressure maintenance during the operation, the puncture guide tube will adopt the way of combining the round-hole silicon sheet and the funnel-shaped slotted-cut silicon sheet, which can respectively bear 4Kpa pressure.

Products - Endo Pouch

15 WickiMed Supply Company LTD

EndoRetrieval pouch

Retriever is placed inside the abdomen during laparoscopic minimally invasive procedures, which can contain organs with certain volume. For easy operation, the retriever is concealed in outer tube after folding and the mouth needs opening after it arrive in the abdominal cavity by puncture outfit guide tube.



Product type	size (inch)	capacity (ml)	length (cm)
Memory type	3*6	100	23
Memory type	3*6	200	23
Memory type	5*7	400	23
Memory type	5*9	600	23
Separation type	3*6	100	29
Separation type	3*6	200	29
Separation type	5*7	400	29
Separation type	5*9	600	29
Simple type	3*6	100	20
Simple type	3*6	200	20
Simple type	5*7	400	20
Simple type	5*9	600	20

Humanized Operation Handle

Ergonomics design makes it able to be operated with one hand and more easy to set multiple retriever concealed design. All joints are put inside the cannula to prevent from scraping or drop out of joints, which are safer and more convenient to operate.



16

Simple Type Retriever

The retriever is protected inside cannula, easier to go through the puncture outfit;

Memory-metal-loop retriever mouth opens automatically, reduce the degree of difficulty;

Convenient to use, cost-effective;

10mm 3*6 and 10mm 5*7

10mm 200ml and 10mm 400ml are optional.



To achieve the effect that the retriever mouth opens automatically, the retriever is fixed with nickel-titanium wire of 0.4mm diameter, convenient for concealed inside the outer tube after folding. In the meantime, the retriever will be brought out when the inside tube go through the outer. Then, the soft and flexible nickel-titanium wire will promptly and automatically open the mouth, convenient for laparoscopic minimally invasive operation.

Products - Suction Irrigation Set

9 WickiMed Supply Company LTD

10

Suction and Irrigation Set

Completely solve the problem of laparoscope smog

To suck and filter the smog completely and immediately is the safest solution!

Operation step 1

Immediately suck and filter new smog;

Use the suction and irrigation set to suck smog during the surgery;

Use the suction function when there is large amount of smog;

Suck and filter the smog;

Use the filtration function when there is a small amount of smog;

◇ High-density activated carbon and double-filtration filter membrane design;

◇ Filter the minimum diameter of 0.03um smoke particles;

◇ No influence on CO2 abdominal pressure when filtering;

◇ Hand adjustment is available when filtering;

Operation step 2

Completely filter abdominal smog;

The puncture cannula connects the filter before the surgery completed;

◇ High-density activated carbon and double-filtration filter membrane design;

◇ Filter the minimum diameter of 0.03um smoke particles;

◇ Hand control is available when filtering;



Products	Cannula DIA (mm)	Steel Tube length (mm)
suction and irrigation set	3、5、10	200、280、330、450



- Self-purification function keeps the suction clear; the large caliber suction can suck blood clot, broken bones slag, muscle and adipose tissue shed off the surgery. No blocking phenomenon!

- Unity of steel tube and main control unit's tube inner diameter;
- Reduce the condition of suction and washing the blocking pipes;
- Improve the washing and suction capacity;
- Anti-skid treatment is easier to be operated;

- Self-locking design enhances security and usability;
- Product structure increases the welding treatment, which improves product sealing;

- The coil spring joining increases the feedback strength and smooth washing and suction. Convenient to use!
- Fast disassembly design helps replacing or adding other equipment more immediate and coordinate;

Expand practicability

The coil spring joining strongly increases the feedback strength and better hand press feeling;
Strong integrity of the product and no leakage phenomenon during the surgery;
Design transparently, simple structure and easily for observation and operation;
More functions are waiting for you...

Products - Suction Irrigation Set

11 WickiMed Supply Company LTD

12

Suction and Irrigation Set

Two kinds of quick-released interfaces are easy to be replaced.



Detachable Interface

Screw-type Interface



Detachable Plug

Screw-type Plug



Replaceable-type Electrosurgery Electrode

This product is available for different surgery operation with different body parts, which is also convenient to be replaced. Safe and high-efficient!



Needle

L-Hook

J-Hook

Spatula

*Optional stealth ER eschar resistant coating



Needle

L-Hook


J-Hook

Spatula

Products - Scissors and Forceps

High-frequency Electrocoagulation Scissors

Ergonomics control handle, accurate operation feeling;
Stainless steel scissors with radians, adaptive to more surgical requirement;
Excellent wire performance, guarantee the surgical fluency;
Rotating degree scissors cannula, more rapid reactivity on complex abdominal environment;



Ergonomics design handle




Stainless steel scissors with radians



33cm-long scissors cannula

High-frequency Electrocoagulation Forceps

Forceps with multiple size, adaptive to different surgical requirement;
High-precision control handle, excellent operating performance;
Rotating degree forceps cannula, more rapid reactivity on complex abdominal environment;



- Extensive range of jaw-tip configurations providing complete surgical offering
- Robust, wide-opening, dual-articulating jaw-tips providing optimal application during a procedure
- A comprehensive range of premium quality, single-use 5mm monopolar instruments with two different handle designs, providing surgeons with the ultimate choice and convenience



Products - OEM

Game-changing Surgical Smoke Pencil
PlumePen® Ultra



Smoke Evacuation Pencil


NISImagine

Making scarless and incisionless surgery possible



Surgical Robotics

Current OEM Products – ESU Pencil



Current OEM Products – Smoke Pencil



WickiMed ODM Products

- Disposable Surgical Mask with Elastic Ear Loop
- 3-ply design (using medical grade melt-blown fabric)
- Comply to EN 14683:2019 Type IIR



Certification

- * ISO 13485:2003
- * ISO 11135
- * CE (Europe)
- * FDA (USA)
- * CFDA (China)
- * PMDA (Japan)
- * MDSAP (Pending)



Regulatory Support

Establishment Registration & Device Listing

[FDA Home](#)
[Medical Devices](#)
[Databases](#)



1 result found for **Establishment Registration**
or FEI Number : 3010601992 Owner Operator
Name : wickimed

[New Search](#)

Establishment Name	Registration Number	Current Registration Yr
Wickimed(Huizhou)Medical Equipment Manufacturing Co.,Ltd. CHINA	3010601992	2020
<ul style="list-style-type: none"> Electrosurgical, Cutting & Coagulation & Accessories 		Contract Manufacturer
<ul style="list-style-type: none"> Electrosurgical, Cutting & Coagulation & Accessories - 3cm TUNGSTEN MICRONEEDLE (Model:138504,138504-10); 4cm TUNGSTEN MICRONEEDLE(Model:138505,138505-10); 5cm 45 Angle TUNGSTEN MICRONEEDLE(Model:138506,138506-10); Electrosurgical Pencil With Electrode,Model:WP200,WP200-B100S,WP200-C100S,WP200-L31S,WP200-N100S,WP200-B90S 		Contract Manufacturer
<ul style="list-style-type: none"> Electrosurgical, Cutting & Coagulation & Accessories - Medline 		Contract Manufacturer
<ul style="list-style-type: none"> Unit, Oral Irrigation - Power Floss Manual Water Flosser (Model: PW011120); Power Floss Manual Water Flosser (Model: PW148); Power Floss Manual Water Flosser P/2 (Model: 29603); Power Floss Manual Water Flosser P/4 (Model: 29616) 		Contract Manufacturer
<ul style="list-style-type: none"> Clamp, Vascular - OMED DISPOSABLE BULLDOG 		Contract Manufacturer
<ul style="list-style-type: none"> Pneumoperitoneum Needle - Veress Needle, Models: WVN0112T, WVN0115T 		Foreign Exporter; Manufacturer
<ul style="list-style-type: none"> Insufflator, Laparoscopic - Veress Needle, Models: WVN0112T, WVN0115T 		Foreign Exporter; Manufacturer
<ul style="list-style-type: none"> Laparoscope, General & Plastic Surgery - Trocar,Model: Auto-Locking Trocar, Bladeless Trocar, Visible Trocar 		Foreign Exporter; Manufacturer
<ul style="list-style-type: none"> Laparoscope, General & Plastic Surgery - Suction Irrigation Set, Model: WSI105280T, WSI105330T, WSI105450T 		Foreign Exporter; Manufacturer
<ul style="list-style-type: none"> Laparoscope, General & Plastic Surgery - Equipment Pouch, Models: WEP010306A, WEP010306B, WEP010507A, WEP010507B 		Foreign Exporter; Manufacturer
<ul style="list-style-type: none"> Instrument, Manual, Surgical, General Use - Cautery Tip Cleaner,Model:118438,DYNJ01208 		Contract Manufacturer
<ul style="list-style-type: none"> Apparatus, Exhaust, Surgical - OR-VAC Electrosurgical Pencil Adapter; SnapEvac Electrosurgical Pencil Adapter 		Contract Manufacturer
<ul style="list-style-type: none"> Electrosurgical, Cutting & Coagulation & Accessories - GoldVac Slim; PlumePen Pro; PlumePen Ultra 		Contract Manufacturer

Facilities & Production Capability

Autoclave Chamber



Projector



Incubation Chamber



Aging Chamber



Flexing Tester



Negative Pressure Tester



Positive Pressure Tester



Sliding Force Tester



Facilities & Production Capability

Class 2 & 3 Production Facilities

10k & 100k Clean Room 1500 sq meter



Facilities & Production Capability

Class 2 & 3 Production Facilities


10k & 100k Clean Room 1500 sq meter



Facilities & Production Capability

Automation

Soldering Process



auto soldering 1.mp4



auto sorting 1.mp4



auto sorting
2.mp4

Stable Quality

Price Competitive

Delivery (Manpower Shortage)

Service and Capabilities

- Laparoscopic Product Production



Packaging and Sealing



Facilities & Production Capability

EO Sterilization Centre (ISO 11135) 15 CBM x4



Facilities & Production Capability



Thank you



medical@wikimed.com
www.wikimed.com